


Accelerating spatial proteomics: Novel on-center TMA core scanning to enable scalable, high-throughput phenotyping for large-cohort TMAs



J Spencer Schwarz¹, Kayla Cashion¹, Charles Jackson¹, Mathew Ingalls¹, Daniel Jimenez Sanchez¹, Joseph Lownik², Oliver Braubach¹

1. Bruker Spatial Biology, St. Louis, MO, USA
2. Cedars-Sinai Medical Center, Los Angeles, CA, USA

Introduction

The etiology of cancer is extraordinarily diverse. Studies that report greater sample variation are more robust, leading to better patient stratification and better treatments. Spatial biology has the potential to improve diagnosis and prognosis of cancer, yet broad adoption of these spatial assays in translational and clinical studies has been hindered by acquisition times and file sizes. The CellScape Precise Spatial Proteomics (PSP) platform generates exceptional data with high signal to noise (SNR) and a spatial sampling of 182 nm/pixel, all while collecting 0.815 mm² of sample area per field of view (FOV). This makes the CellScape platform an attractive tool that balances data quality with acquisition speed and overall instrument utility. Yet, there are cases in which a project may require still faster data collection.

Here we demonstrate the CellScape platform's FalconFAST mode, a primary 10x 0.45NA objective (10x 0.45NA WD 4.0mm Plan Apo Lambda D, Nikon Instruments), to increase data collection speed. In this configuration, the instrument can record 3.18 mm² per FOV and quickly scan through large biopsies and tissue microarrays (TMAs). We qualified these data by demonstrating fulfillment of three requirements fundamental to any multiplexed immunofluorescence (mIF) data pipeline. The data must be suitable for 1. qualitative marker scoring; 2. automated image segmentation; and 3. automated phenotyping.

Method

The CellScape workflow centers around EpiclF™ signal removal technology. EpiclF technology effectively and rapidly deactivates fluorophores while leaving tissues and antibodies unaffected. Using this workflow with the modified CellScape platform, we deployed VistaPlex™ Spatial Immune Profiling and Tissue Architecture Assay Kits on FFPE tissues. The 34 pre-validated assay markers provided the information needed for cell segmentation, immune cell identification and architectural analysis of the stroma, epithelia, and the vasculature.




Figure 1. Experimental workflow using CellScape Precise Spatial Proteomics and EpiclF technology. Cycles of staining, imaging, and signal removal detect biomarkers with spatial context at single-cell resolution. Signal removal facilitated by filtered photobleaching and the EpiclF Solution to provide a safe, gentle, and effective fluorophore removal.

Scan here to download or learn more

Figure 2. Marker qualification and staining evaluation. Panel a, and inset, show major histological features of the tonsil: tonsillar crypt, 1, germinal center, 2, and interfollicular zone, 3. Panel b shows cell segmentation markers with sufficient resolution and contrast to achieve accurate segmentation (see also Fig. 4). Panels c-g illustrate epithelial boundaries, c, cell-cell interactions, d, distinct but closely packed stromal, e, and immune cells, f, and sub-cellular features, g. h shows magnified views of all markers present in the dataset. Right panels and inset show markers collected with the standard 20x 0.80NA objective (tissue and assay matched), while the left panel and inset show data collected with CellScape FalconFAST mode.





Figure 3. Field-of-View Placement comparison. FOV layouts for FalconFAST scanning, a, and the CellScape platform's standard scanning mode, b. Data in c show the time required to collect each data set, and the resulting file size of the data. The 0.305 FOV fraction is reflected in the per-marker scanning estimate. A pixel fraction of appx. 1/4th is expected from the difference in sampling rates. Data size (pixel count) impacts image processing time.



Results

Figure 4. Single-marker and VistaPlex Assay Kit qualitative scoring. Each marker was manually scored to evaluate marker specificity and sensitivity. Data below show that most markers pass the > +1.5 threshold that is required in our quality control. Scores are applied to each stain on a scale from -2, non-specific staining, through 0, no staining, to +2, appropriate staining. The radar plot below represents the average of two independent assessments.

